# multisportpro

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

## MULTISPORTPRO (4 CHANNELS) DESCRIPTION OF THE UNIT

- 1 Left key : Back or select another choice
- 2 OK key : Validation after each selection
- 3 C key (clear) : Back after validating
- 4 Right key : Choice selection key (programs, muscular zones.)
- 5 ON/OFF key : Turns the unit On and Off / emergency stop
- 6 4 Power key :  
+ : Increase power / level 1 to 30  
- : Decrease power / level 30 to 1

## DESCRIPTION OF THE SCREEN

- Battery level
- Pg 88** Displays the program number chosen
- Displays the countdown
- ALT** Alternate muscles option (choice between programs PG01 to PG08). Train 2 muscle groups alternately.
- SIM** Simultaneous muscles option (default choice). Train 2 muscle groups simultaneously.
- Enable/disable warm-up phase
- Sport program
- Health program

0459 (directive 93/42/CEE)

Made in P.R.C.

## PROGRAMS

### SPORT PROGRAMS

#### SKIN PREPARATION

Lowers the impedance of the skin barrier. Prepares your skin for the currents emitted by the programs.  
**Duration :** 2 min active phases

#### RESTORE MUSCLE VOLUME

This program prevents and fights loss of muscle volume resulting usually from a lack of use of a limb following a period of immobilisation.  
**Duration :** 40 min active phase + 6 min recovery

#### REINFORCE MUSCLE VOLUME

Apply this program only on muscles with a normal volume. PG04 and PG05 are used to develop muscle strength.  
**Duration :** 20 min active phase + 6 min recovery

#### MUSCLE RESISTANCE

The stimulation frequencies aim at increasing the ability to sustain major muscle efforts.  
**Duration :** 30 min active phase + 6 min recovery

#### MUSCLE ENDURANCE

After using PG04 to acquire the ability to sustain major efforts, use program PG05 to extend high intensity muscle efforts over time by maximum stimulation of slow fibres.  
**Duration :** 30 min active phase + 6 min recovery

#### MUSCLE FITNESS

Use this program to maintain the level of muscle endurance, resistance, and volume acquired.  
**Duration :** 30 min active phase + 6 min recovery

#### EXPLOSIVE FORCE

This program is a logical extension of PG03 and PG04. It allows reaching maximum muscle force by stimulating very high speed fibres, with maximum contraction intensity.  
**Duration :** 30 min active phase + 6 min recovery

#### RECOVERY

Using a frequency range between 0.5 and 8 Hz, this program ensures muscles recover faster through relaxing, endorphin-based effects, and increases blood flow.  
**Duration :** 20 min active phase

### HEALTH PROGRAMS

#### MASSAGE

This program uses frequencies modulated between 3 and 10 Hz for optimum massaging effect.  
**Duration :** 20 min active phase.

#### MUSCLE RELAXING

The frequency range between 0.25 and 3.5 Hz of this program accelerates muscle relaxing.  
**Duration :** 20 min active phase

#### ACTIVE RECOVERY

Use this program immediately after performing intense efforts. All frequencies generated by PG11 maintain muscle work after effort, enhancing blood flow and the release of endorphins. Each of these frequency sequences is alternated with a relaxing effect.  
**Duration :** 20 min active phase

#### TENS

This program operates by cross-skin nerve stimulation, delivering a range of pulse frequencies in the local area where the electrodes are placed.  
**Duration :** 20 min active phase

#### ENDORPHINIC

This program relieves muscle tensions effectively.  
**Duration :** 20 min active phase

#### VASCULARIZATION

This program relieves numb legs sensation.  
**Duration :** 26 min active phase

## BEFORE USING THIS UNIT

**You have just acquired an effective and powerful professional unit.**

**Take your time to discover how to use it.**

- Read the instructions carefully
- Before using it, charge the unit (photo I. E page 3)
- Test your unit on your quadriceps using the quick test procedure below

## OPERATION QUICK TEST

### TESTING YOUR UNIT ON YOUR QUADRICEPS

- Take a ½ belt, apply cream on the built-in electrodes (photo III. E page 4).
- Connect it to the first channel (match the colours of the connectors, 2 by 2).
- Start the unit using the  key and press the "OK" button 6 times.
- Increase power on this channel gradually.
- This discovery procedure shows your capacity in using the unit's level of performance..

### \* THERAPEUTIC APPLICATIONS

- **Private:** to strengthen and / or to maintain muscle strength.
- **Physiotherapists:** to ease muscular tensions.
- **In hospital environment:** to prevent muscular atrophy and to treat pain.

## I. CHARGING THE BATTERIES

The unit contains rechargeable batteries. Before using it, charge it during 12 hours by connecting it to mains. Do not use this unit while batteries are charging. When it is not used, charge the batteries for 12 hours every 4 months. When the batteries are empty, the unit switches to standby automatically, protecting them against untimely discharging. To charge the unit, connect the charger supplied with the unit (photo I. E page 3)

### BATTERY CHARGE DISPLAY

When the battery is charging, the battery indicator displays , , then  to show the current charging status. The charging time countdown starts, from  to . Once charged, the number  appears in the left angle of the timer pictogram. The battery pictogram displays . When the charging is complete, unplug the adapter from the electrical outlet to insulate electric source. When the charger is disconnected, the screen switches off. When the battery indicator is empty , this means that it is discharged.

## II. STARTING THE UNIT

### STEP 1 : DISPLAY - POWER ON

Beware while starting the device ( red button  )  
Push briefly (maximum 1 second ) on .  
All pictograms appear and blink for 3 seconds, then the front muscles blink.  
If you make a long push, your screen flashes and indicates you the totality of the pictograms which it contains.  
To stop the device, push briefly on .  
Or to return in the menu made 2 brief pushes on .

### STEP 2 : SELECTING THE MUSCLE ZONE(S) TO STIMULATE

#### Choosing the 1<sup>st</sup> zone :

The screen displays , meaning you can choose a first muscle zone.  
- To choose the front muscles, validate with the "OK" key.  
- To choose the back muscles, press the  key.  
- To return to the front muscles, press the  key.  
If you have chosen the front muscles : the 1<sup>st</sup> muscle zone at the bottom blinks: quadriceps.  
If you have chosen the back muscles : the 1<sup>st</sup> muscle zone at the bottom blinks: gastrocnemius (calves)  
Press the  key to move upwards and scroll the different muscle zones.  
To go back, press the  key.  
After selecting the muscle group, validate using th "OK" key

### Choosing the 2<sup>nd</sup> zone :

All pictograms blink and the screen displays 02.

The screen proposes to choose a 2<sup>nd</sup> muscle zone :

To choose a 2<sup>nd</sup> muscle zone, validate again either the front muscles or the back muscles, validate "OK".

Choose again another muscle zone, then validate "OK".

If you wish one muscle zone only, press the "◀" key and go to step 3.

The muscle zone(s) selected is (are) displayed.

### STEP 3 : SELECTING THE PROGRAM

Programs on the left (SPORT programs) blink.

If you choose the SPORT programs, validate "OK".

To choose the HEALTH programs, press the "▶" key and validate "OK".

Scroll through the programs by pressing the "▶" key.

The screen first displays the 1<sup>st</sup> lower pictogram while displaying the corresponding program number at the top.

To go back, press the "◀" key.

You have now chosen the muscle zone(s) to be stimulated, and the program.

### STEP 4 : CHOOSE THE "WARM-UP PHASE" AND ALT SIM OPTION

(Only for programs P602 to P607).

The screen proposes a warming-up phase before starting the active phase of the program.

The "WARM-UP PHASE" pictogram blinks in the top right corner.

To select it, press "OK", otherwise press the "◀" key, and validate "OK" (it disappears from the screen).

Then the "SIM" pictogram blinks.

To select "ALT", press the "◀" key, to select "SIM", press the "▶" key, then validate "OK".

If you have selected the warm-up phase, "SIM" is displayed when starting the program, as the warm-up phase is only in the simultaneous mode.

The screen is displayed :

- the battery level,

- phase duration (countdown in minutes),

- timer pictogram,

- the "ALT" or "SIM" option and warm-up phase if you have selected them.

The pictogram for the program chosen is lit, the muscle zone(s) chosen blink alternately.

### STEP 5 : LAUNCHING THE START - ADJUSTING INTENSITY

Then, press the "+" buttons of the intensity keys matching the wires connected to your electrodes or your belt to start the program.

If the battery is low, do not start the program.

The power level(s) are displayed at the bottom of the screen. The intensities vary during the session according to the program's content.

A long press on the "—" key returns the power rapidly to 0.

### STEP 6 : PROGRAM PERFORMANCE :

**Use of the "C" (Clear) key :**

During the program, press the "C" key during a phase to set its time to zero.

The next phase starts unless you have reached the end of the program.

In this case, this key allows returning to selecting the 1<sup>st</sup> muscle zone (as at normal end of program).

### Disconnecting the electrodes during the program :

If an electrode is disconnected, the channel concerned is reset to 00 and blinks for 30 sec, before returning to a fixed 00.

During these 30 sec, the electrode(s) can be reconnected, and power restored to the level required.

### STEP 7 : END OF PROGRAM - STOPPING

At the end of the session, the screen proposes to start a new session by blinking the front muscles to select the 1<sup>st</sup> muscle zone.

Restart the starting procedure if you wish to continue on another muscle zone.

Otherwise, press the "⏻" key to stop.

The unit switches to standby after 5 minutes idle.

### Recommendations :

The pictograms for upper muscles and for lower muscles are different. Avoid selecting upper muscles and lower muscles, as the unit would give priority to upper muscles, reducing somewhat the yield of work on lower muscles.

## III. USING THE UNIT

### A. ELECTRODES FUNCTION

This function allows using the adhesive electrodes by connecting them to the unit

#### STEP 1 :

Take the electrodes and remove carefully their protection, dampen them one by one with water on your fingertips on the adhesive gel, and wait for 2 minutes (diagram I. A, page 3).

#### STEP 2 :

Fit the electrodes on your body as shown in diagram II. A or B or C, page 3 : Positions are shown for information only

**Important :**

For optimal use, carefully read each page and look at the examples of use (photos II page 3).

Examples of use (photos II page 3)		
B = connector channel B (middle blue) ; C = connector channel C (light blue) ; SUP = upper limbs ; INF = lower limbs		
A photo - Chest	B photo - Shoulders	C photo - Trapezius
1 - Large chest ⚠ Do not place the electrodes close to the heart	2 - Medium deltoid 3 - Posterior 4 - Anterior	9 - Upper trapezius 10 - Medium trapezius

#### STEP 3 :

Take connector (CEV8) (photo I. B, page 3) and connect it :

1° - to the electrodes (photo I. C, page 3)

2° - then to the unit (photo I. D, page 3).

#### STEP 4 :

The unit is connected to the electrodes through the CEV8 connector, made of 2 groups of 4x 1m. long wires terminated by 2 groups of 2 connecting terminals.

Fasten your MULTISPORTPRO to the back of your clothing or underwear using its clip.

#### STEP 5 :

Start your unit following the "Starting the unit" procedure described in paragraph II.

Contractions should be firm but never painful. If you feel an unpleasant sensation around the electrodes : bring the power-on key back to minimum, move the electrode(s) and gradually increase the power.

Once the adjustments have been optimised in this way, leave your appliance to work for the full duration of the program.

#### STEP 6 :

Once the program is complete, disconnect the electrodes from the (CEV8) connector, then the latter from the unit.

Then pull off the electrodes gently, dampen them slightly with a few drops of water, then return them carefully onto their support.

Store the adhesive electrodes in a fresh location, in a plastic bag (e.g. vegetable container in the refrigerator).

### B / BELT FUNCTION

The multifunction belt can be used complete, or separated into 2 to form 2 half-belts.

The full belt must be used on the stomach, buttocks, lower back areas.

Once unfastened, the 2 half-belts can be used on the arms, thighs, or calf.

The multifunction belt has built-in electrodes on which the SPORT-ELEC® cream is applied.

This cream enhances contact between the electrodes and the skin.

Its massaging and electro conducting properties improve comfort of use and enhance stimulation.

The cream can be also used by massaging on the body.

If you want to buy another SPORT-ELEC® cream, please order directly by going to our online shop : [www.sport-elec.com](http://www.sport-elec.com), or contact your local supplier.

**You may not use the belt without the cream provided.**

**Do not use any other cream than that the one provided by SPORT-ELEC® Institut.**

#### SAFETY PRECAUTIONS BY USING THE CREAM :

- In case of allergy, stop to apply the cream.
- Do not swallow.
- Do not leave with the range children.
- In case of contact with eyes rinse with water.

#### PRESENTATION OF THE BELT :

##### • Front side (photo III. A, page 4) :

- 1 - Central fasteners    2 - Wire storage    3 - Connection terminals    4 - Scratch strip

##### • Electrodes side (photo III. B, page 4) :

- 4 built-in electrodes

#### STEP 1 : CONNECTING AND FITTING THE BELT.

- Check your MULTISPORTPRO unit is STOPPED
  - Take the (CEV8) connector (photo I. B, page 3) and connect it :
    - To the belt connectors per groups of colour (photo III. C page 4).
    - Then to the unit ( photo III. D, page 4).
  - Apply a thick layer of electro conductive massage cream onto the electrodes built into the belt (domed part) (photo III E page 4)
  - Place the belt or 2 half-belts onto the zone(s) to be stimulated, and adjust them using "scratch" strips to ensure the built-in electrodes are properly in contact with the skin  
Photos IV page 4 : A - stomach, B - back, C - arms, D - thighs, E - calf, F - back.
  - Move the belt slightly on left or right (with device powered off) to feel the contractions on the abdominal area.
- If an electrode hurts or if you feel any tingling, stop the device, and check the electrode is properly smeared with cream.

##### TIP :

- Fitting the full belt :  
To well-positioned the full belt at best, cross the 2 strips to flip the 2 ends on the outer scratch part (photos V. A page 4).
- Fitting the half-belt on the arm :  
Adjust at the belt mid-size, according to the outline of your arm before positioning it, then start over if necessary to ensure proper adjustment (photos V. B page 4).

#### STEP 2 : STARTING THE BELT FUNCTION

Once your unit is connected, install yourself comfortably and start your unit using the procedure "Starting the unit" detailed previously in paragraph II.

#### STEP 3 : STOPPING THE BELT FUNCTION

Once the program is completed, disconnect the unit from the (CEV8) connector, and the latter from the belt. Remove the excess cream on the built-in electrodes using a wipe dampened with water, then dry them.

## IV. RECOMMENDATIONS OF USE

The unit is recommended for physiotherapists for functional re-education.

- You can restart a program's cycle on another muscle zone or change the program while remaining on the same muscle zone, without exceeding however 40 min use per muscle zone and per day.
- For optimum use, refer to the diagrams shown in this leaflet.
- Clean the unit and its cables with a cloth dampened with a disinfecting solution (Mercryl type) diluted in water at 1/10<sup>th</sup>.
- The adhesive electrodes can be rehydrated before and after use with a disinfecting agent (e.g. Hexamedine).
- The adhesive electrodes must be repositioned on their support, and stored in a fresh storage between +5°C and +10°C.
- The adhesive electrodes have a lifetime of about 40 uses, you can order them directly from the online shop [www.sport-elec.com](http://www.sport-elec.com), or from your reseller.

#### PRECAUTIONS OF USE :

- Do not use the unit on your face
- Do not leave within reach of children

##### Ask for advice from your physician :

- if you have electronic or metallic implants (pins, etc.)
- in case of skin condition (wounds,...)

#### WARNINGS :

- Keep away from children or people with mental disabilities.  
The device contains small pieces that can be swallowed.
- In order to prevent any risk of strangulation, please keep the electrodes wires away from children.
- If the device is not likely to be used for some time, the batteries should be removed
- Multisportpro is meant for personal use. For reasons of hygiene, the electrode should only be used on one person.
- Do not use the device by driving or by working on a machine.
- Do not use this electro stimulation unit in a damp environment, a bathroom, a sauna, or close to electric machinery connected to mains and earthed, or even piping earthed.
- Use this device in a clean place (without dust nor dirt...).
- The simultaneous connection of a patient to a high frequency surgical unit may cause burns at the points of contact of the stimulator's electrodes, and the stimulator may be damaged.
- Operation in close proximity of a short wave therapy unit may cause instability in the stimulator's output power.
- Do not, under any circumstances, attempt to open or repair the Multisportpro yourself! Only have repairs carried out by the customer service department.

#### TROUBLE SHOOTING

If your device doesn't work properly, do not use it and contact your retailer.

#### RECOMMENDATIONS :

-  Please apply regulations enforced in your country when disposing your product at its end of life.
-  Therefore, please discard it in a location planned to this effect, ensuring it is disposed in an environment-friendly manner.

#### ⚠ CONTRA-INDICATIONS : Do not use SPORT-ELEC® :

- on the heart area
- If you carry a pacemaker
- If you are pregnant
- during first 6 weeks after baby birth
- In case of neuromuscular diseases
- In case of hemorrhagic diseases
- On the face and skull

#### RECOMMENDATIONS FOR MAINTAINING AND CLEANING THE BELT :

IMPORTANT : NEVER wash in machine and under running water. The dryer and ironing facilities are forbidden.

- Do not dry-clean
- Do not expose to bleach
- Clean the belt using a cloth dampened with a detergent, then wipe with a dry cloth.
- Clean the unit with a tissue dampened with a disinfecting solution (Mercryl type) diluted with water at 1/10<sup>th</sup>.
- Clean the electrodes with a mix of water/hypoallergenic detergent, then rinse thoroughly

**CUSTOMER SUPPORT: +33 (0)2 32 96 50 50**

## V. TECHNICAL SPECIFICATIONS

### UNIT :

Class IIa medical device in accordance with Directive 93/42/CEE	
Software version.....	E
Type of current.....	Dual phase
Number of programs.....	8 "SPORT" programs - 6 "HEALTH" programs
Number of program variations (or phases).....	344
Number of adjustable independent outputs.....	4
Maximum intensity on a 500 $\Omega$ load/channel.....	100 mA / 500 $\Omega$ / channel
Frequency range.....	0,25 to 160 Hz
Pulse width range.....	70 to 360 $\mu$ s
Power supply.....	4.8 V - 1300 mAh - NI-MH Rechargeable battery
Mains supply.....	THA239025 : 230VAC 50Hz 100mA - 9VAC 250mA RSS1002-055110-W2B-B : 100-240VAC 50/60Hz 0.18A - 12VDC 250mA
	The adapter is part of the equipment. Only use mains power after reticulation by the power unit.
Unit dimensions and weight.....	80 x 160 x 35 mm - 340 gr
Automatic stop at end of each program.....	Yes
Activation safety.....	100%
Temperature for stocking / Relative humidity.....	0°C to + 45°C / 10% to 90%
Temperature for use / Relative humidity.....	5-40°C / 20-65%
Caution !.....	 Read the manual thoroughly
Guarantee of the unit / Lifecycle.....	2 years

### IP22

Device in contact with the patient, excluding the heart muscles.

Protected against the ingress of foreign object  $\geq 12.5$  mm & against dripping water when tilted up to 15°.

### RoHS

These units are free of materials likely to have a negative impact on the environment.

### BELT :

Size.....	Universal size (about 150 cm long) (i.e. 2 half-belts 75 cm)
Weight.....	150 g
Fabric.....	Polyester laminated on foam
Built-in electrodes.....	Carbonated silicon
Estimated lifetime.....	around 10 months
Class I medical device in accordance with Directive 93/42/CEE	

### ELECTRO CONDUCTIVE CREAM :

Packaging.....	75 ml / 2.54 fl Oz
Manufacturer reference.....	VC 57B/53-148
Class I medical device in accordance with Directive 93/42/CEE	

### ELECTRODES :

Electrodes.....	Adhesive, hydrophilic and hypoallergenic
Size.....	89 x 50 mm
Estimated lifetime.....	40 uses
Class I medical device in accordance with Directive 93/42/CEE	

### ACCESSORIES GUARANTEE :

The belt, cream, and electrodes are designed for personal use. For reasons of hygiene, there will be no trade-in nor exchange.

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### EN - WARRANTY

This appliance is guaranteed for 24 months from the date of purchase. During the period of the guarantee, the distributor will take responsibility for defects occur, the appliance will be repaired or exchanged as the distributor deems necessary. Damage arising from improper use or normal wear and defects, which do not materially affect the performance of the appliance, is not covered under this guarantee. The guarantee is not valid if unauthorized repairs have been attempted or if replacement parts used are not approved by the supplier. The guarantee is valid only if it is stamped and signed by the retailer at the time of the purchase. For all claims made during the period of the guarantee, send the appliance and the warranty card to your retailer. All other claims are excluded from the guarantee unless specifically prohibited by law.

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